

®

3V AGILE

Central Venous Catheter Kit



SEV®
VASCULAR TECHNOLOGIES

®

3V AGILE

Central Venous Catheter Kit Accessories



Syringe




Guide wire



Injection Cap (S)



Dilator



Fixator wing



Scalpel



Introducer Needle

SEV®
VASCULAR TECHNOLOGIES



3V AGILE

Central Venous Catheter Kit

Lumen Type	Available Sizes		
Single Lumen	16g/16cm	18g/8cm	20g/6cm
Double Lumen	7Fr/16cm	4Fr/6cm	8Fr/8cm
Triple Lumen	7Fr/16cm	4.5Fr/6cm	5.5Fr/6cm
Quad Lumen	8.5Fr/16cm		

Contact Information:

S3V Vascular Technologies Limited
Plot No. 65-D4, Hootagalli Industrial Area,
Hootagalli, Mysore- 570018, Karnataka, India.
Tel: +91 8212410016
Email: info@s3vvascular.com
Web: www.s3vvascular.com

